# Heraeus

## **Heramat**

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<b>14</b>	<b>Dokumentenhistorie.....</b>	<b>D 34</b>
<b>Operating Instructions - Heramat C</b>		
	<b>English .....</b>	<b>GB 35</b>
<b>1</b>	<b>Scope of applicability .....</b>	<b>GB 35</b>
1.1	General.....	GB 35
1.2	Machine name and type .....	GB 35
1.3	EC - Declaration of Conformity .....	GB 35
1.4	Relevant EC directives .....	GB 35
<b>2</b>	<b>Instructions on using the unit safely.....</b>	<b>GB 36</b>
2.1	Explanation of symbols.....	GB 36
2.2	Transport damage .....	GB 36
2.3	Operator's obligations.....	GB 36
2.4	Unit diary / operating manual.....	GB 36
2.5	Safety instructions .....	GB 37
<b>3</b>	<b>Intended use .....</b>	<b>GB 37</b>
3.1	Working rules.....	GB 38
<b>4</b>	<b>Scope of delivery.....</b>	<b>GB 38</b>
<b>5</b>	<b>Structure and function.....</b>	<b>GB 39</b>
5.1	Heramat C - overall view .....	GB 39
5.1.1	What do the LEDs in Fig. 1 Overall view of the Heramat C mean?.....	GB 39
5.2	View of control panel .....	GB 39
5.2.1	What do the LEDs in Fig. 2 View of operating panel mean? .....	GB 40
5.3	Rear view.....	GB 40
5.3.1	What do the LEDs in Fig. 3 Rear view of machine mean? .....	GB 40
5.4	View of combustion table.....	GB 41
5.4.1	What do the LEDs in Fig. 4 View of combustion table mean?.....	GB 41
5.5	View of accessories .....	GB 41
5.5.1	What do the LEDs in Fig. 4 View of combustion table mean?.....	GB 41
<b>6</b>	<b>Unit description .....</b>	<b>GB 42</b>
6.1	Description of keys, fields and displays .....	GB 42
6.1.1	Data keys.....	GB 42
6.1.2	Control keys.....	GB 42
6.1.3	Program keys and fields .....	GB 42
6.1.4	Displays .....	GB 42
6.2	Description of acoustic signals .....	GB 43
<b>7</b>	<b>Putting the unit into operation .....</b>	<b>GB 43</b>
7.1	Transport .....	GB 43
7.2	Installation .....	GB 43
7.2.1	Vacuum pump connection .....	GB 43
7.2.2	Furnace connection to the mains power supply .....	GB 43

7.3	Initial operation .....	GB 44
7.3.1	Initial combustion process of the furnace .....	GB 44
7.4	Operating the machine .....	GB 44
7.4.1	Switching on .....	GB 44
7.4.2	Program mode.....	GB 44
7.4.3	Set-up option .....	GB 44
7.4.4	Idling mode .....	GB 45
7.4.5	Program selection.....	GB 45
7.4.6	Programming mode .....	GB 46
7.4.7	Entering and altering parameter data .....	GB 46
7.4.8	Description of parameters and their limit ranges .....	GB 46
<b>8</b>	<b>The combustion process .....</b>	<b>GB 49</b>
8.1	Starting the combustion process .....	GB 49
8.2	Idling mode .....	GB 50
8.3	Interrupting the combustion process .....	GB 50
<b>9</b>	<b>Additional features .....</b>	<b>GB 50</b>
9.1	Time remaining to end of combustion process .....	GB 50
9.2	Post Fire .....	GB 50
9.3	The HOLD key .....	GB 50
9.4	Altering parameters during the combustion process .....	GB 50
9.5	Monitoring the vacuum level .....	GB 51
9.6	Night mode .....	GB 51
9.7	LIFT key .....	GB 51
9.8	STOP key .....	GB 51
9.9	Temper .....	GB 51
9.10	Diagnosis mode .....	GB 52
<b>10</b>	<b>Servicing .....</b>	<b>GB 52</b>
10.1	Tests .....	GB 53
10.2	Preheating .....	GB 54
10.3	Vacuum calibration .....	GB 54
<b>11</b>	<b>Technical data .....</b>	<b>GB 57</b>
11.1	General.....	GB 57
11.2	Transportation and storing conditions: .....	GB 57
11.3	Operating conditions:.....	GB 57
11.4	Rating plate .....	GB 58
11.5	Disposal of old equipment according to WEEE .....	GB 58
<b>12</b>	<b>Circuit diagrams .....</b>	<b>GB 59</b>
12.1	Terminal allocation scheme .....	GB 59
<b>13</b>	<b>Heraeus Kulzer adresses worldwide .....</b>	<b>GB 60</b>
<b>14</b>	<b>Document history .....</b>	<b>GB 61</b>

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#### 1 Scope of applicability

##### 1.1 General

Combilabor® is a registered trademark of Heraeus Kulzer GmbH®  
 Author 42415SCH

**These operating instructions apply for:**

Order no.	Type and equipment	Edition
660 02 330	Heramat C, with new CPU from software V 2.00	04.09.2006 15:27

##### 1.2 Machine name and type

Machine name:	Machine type:	Machine no.:
Heramat C	Combilabor CL Heramat C	202 01 001 ff.

##### 1.3 EC - Declaration of Conformity

We, Heraeus Kulzer GmbH, Grüner Weg 11, 63450 Hanau, hereby declare that the machine named below is designed and built to meet the basic relevant health and safety requirements specified in the EC Directive, as is the model launched by us.

If the machine is changed without our approval, this declaration will cease to be valid.

##### 1.4 Relevant EC directives

Machine Directive 89/392/EWG, Appendix II A

EG Low Voltage Directive 73/23/EWG




Electromagnetic Compatibility 89/336/EEC

The unit meets the following safety requirements:

- EN 61010-1 Edition 8/2002
- EN 61326-1

## 2 Instructions on using the unit safely

### 2.1 Explanation of symbols

Symbols	Accompanying word(s)	Explanation
	<b>CAUTION!</b>	<i>Safety-related chapter and sections within the operating manual</i>
	<b>NOTE!</b>	<i>Notes within the operating manual on how best to use the unit</i>
	<b>HOT SURFACE!</b>	<i>Hot surface Danger of burning</i>
	<b>ELECTRIC SHOCK!</b>	<i>Caution electric shock Risk of fatal injury if instructions are not followed.</i>
	<b>DISCONNECT FROM MAINS!</b>	<i>Danger of electric shock if unit is opened Unplug the unit before opening it.</i>
	<b>CHANGE!</b>	<i>Important: Changes have been made to this paragraph. Please read carefully.</i>

### 2.2 Transport damage

Check the machine for damage after receiving it. If it is damaged, report this to the transportation company within 24 hours of delivery. Under no circumstances, work with a damaged unit.

### 2.3 Operator's obligations

In addition to complying with the statutory regulations specified by the manufacturer, the operator must ensure the statutory obligations are observed and implemented in the workplace, i.e. he must train his personnel and comply with industrial safety legislation and any other regulations or laws in force.

For working on and with the unit, the operator must draw up written instructions in understandable form and give these to his employees in their own language. These instructions must be based on the operating manual and written in light of the work to be performed.

### 2.4 Unit diary / operating manual


We recommend that you keep a "machine diary" and document all tests, and major work (e.g. servicing, changes) in it. Keep this operating manual in a safe place so that you can consult it whenever you need safety instructions and important user information.

## 2.5 Safety instructions

The safety instructions for personal safety, the surroundings, and the goods being machined depend heavily on the behaviour of the people using the machine.

Before putting the unit into operation, read the operating manual carefully and follow the instructions to avoid errors and damage, especially damage to your health.

In addition to the instructions in this operating manual, comply with the national laws, regulations, and directives in your country when setting up and operating the unit.

***ELECTRIC SHOCK!***

***Before using the unit, check the mains lead and plug for damage. If they are damaged, do not connect the unit to the mains.***

***HOT SURFACE! The metal surface around the combustion chamber and in particular under the combustion chamber becomes hot during operation, which results in the risk of burning. Avoid direct contact with bare fingers on this surface.***

***Use the tweezers provided.***

***Exhaust gases generated during work must be discharged safely to the atmosphere. The relevant national environmental protection regulations must be observed.***

Work on the unit's electrical equipment must only be performed by trained **service personnel**. The unit must be switched off when such work is performed.

Only use authorised original spare parts and accessories. Avoid other parts as they hide unknown risks.

The functional efficiency and safety of the unit is only guaranteed if the necessary tests, service, and repair work is performed by the **Heraeus service partner** or by personnel trained by us.

Heraeus Kulzer GmbH will accept **no liability** if the unit is damaged by faulty repair work which was not performed by the Heraeus service partner or by trained service personnel, or if original spare parts and accessories are not used to replace parts.

## 3 Intended use

The Heramat C ceramic furnace is a piece of laboratory equipment for producing ceramic veneers in dental technology.

We recommend the use of original Heraeus-Kulzer alloys and materials of the "Heraceram" series. If materials from other companies are used, the user must ensure that they correspond to the state of the art and will have no negative effect on material and equipment.

If, in the event of damage, a causal relationship between other materials used and the technical failure/damage/complaint is determined, a warranty claim against Heraeus Kulzer will not be accepted.

The unit is suitable for setup and operation in the following areas: laboratories, e.g. in business, industry, schools, universities, hospitals etc.

The unit is intended for continuous operation. It is designed for a service life of 10 years.

The unit must be operated by dental-laboratory technicians or assistants who:

- have acquired their skills through technical training or instruction,
- are qualified to perform their duties correctly,
- can detect possible dangers when working with the unit.

### 3.1 Working rules

**CAUTION**

*The unit must not be used for the following purposes.*

- Do not use the unit to warm food.
- Do not use the unit for work during which combustible gases and vapours are released, or for work involving burning with air or in which a dangerous explosive mixture could be formed.
- The unit is also not suitable for the heat treatment of dangerous dusts and fibres.
- Do not operate the unit in non-ventilated rooms.

**CAUTION**

*If particular problems not treated in enough detail in this operating manual arise, please contact the supplier to ensure your own safety.*

**CAUTION!** *Remove any jewellery and put on the necessary personal safety clothing, e.g. hand, face, and body protective clothing.*

## 4 Scope of delivery

Check that all components are in perfect condition on delivery of the unit. If you wish to make a complaint, contact your supplier.

Delivery unit:

- Heramat C unit
- Power cord
- Vacuum hose
- Combustion table
- Tweezers
- Fuses
- Wrench
- Operating manual
- Thermo-tray with 8 pins

## 5 Structure and function

### 5.1 Heramat C - overall view


Fig. 1 Overall view of the Heramat C

#### 5.1.1 What do the LEDs in Fig. 1 Overall view of the Heramat C mean?

- 1 Heramat C unit
- 2 Control
- 3 Combustion table

### 5.2 View of control panel


Fig. 2 View of operating panel

### 5.2.1 What do the LEDs in Fig. 2 View of operating panel mean?

- 4 On/Off switch
- 5 Start/Stop switch
- 6 Program step up/down key
- 7 DATA Display
- 8 PROG (program selection)
- 9 HOLD (program hold)
- 10 LIFT (Lift platform movement)
- 11 ENTER
- 12 CLEAR (delete)
- 13 P:FIRE (afterburner)
- 14 VAC (vacuum query)
- 15 TIME (Time remaining to end of combustion process)
- 16 10-digit keypad
- 17 Parameter field
- 18 Program display

### 5.3 Rear view



Fig. 3 Rear view of machine

#### 5.3.1 What do the LEDs in Fig. 3 Rear view of machine mean?

- 20 Master switch (illuminated)
- 21 Opening for wrench
- 22 Vacuum pump fuse
- 23 Muffle securing element
- 24 Grommet for vacuum hose
- 25 Vacuum pump mains connection
- 26 Ceramic furnace mains connection
- 27 Calibration input
- 28 RS 232 interface
- 29 Rating plate


## 5.4 View of combustion table

### 5.4.1 What do the LEDs in Fig. 4 View of combustion table mean?

- 32 Combustion table
- 33 Lift platform
- 34 Tray

## 5.5 View of accessories

### 5.5.1 What do the LEDs in Fig. 4 View of combustion table mean?

- 35 Power cord
- 36 Thermo-tray
- 37 Object pins
- 38 Spare fuses
- 39 Lift wrench
- 40 Tweezers

## 6 Unit description

### 6.1 Description of keys, fields and displays

#### 6.1.1 Data keys

Key	Description
NUMBERS	Numbers 0-9
ENTER	Confirms selection of a program or parameter.
CLEAR	Deletes the current entry

#### 6.1.2 Control keys

Key	Description
ON / OFF	Switches furnace ON/OFF. Switches to idling mode.
LIFT	Raises and lowers lift platform.
START	Starts a combustion cycle.
STOP	Interrupts a combustion cycle.
HOLD	Temporary "freezing" of a combustion cycle or lift movement.

#### 6.1.3 Program keys and fields

Keys/fields	Description
PROG	Sets furnace to program selection mode
STEPS keys	Allows operator access to desired parameter for programming.
Parameter field	13 parameters that may form the programming of the combustion process are accessible via the [STEP] key (up/down key).
VAC	The word VAC is displayed in the lower section of the parameter field if there is a VACUUM in the combustion chamber.
END	The word END is displayed at the bottom right of the parameter field at the end of the combustion process.

#### 6.1.4 Displays

Displayed data may represent temperature, time or vacuum values.

The program number indicates either the selected program 0-99 or the space / delete symbol [ \_ \_ ].

The symbol [°C] is displayed together with the temperature if the user uses the Celsius unit; Fahrenheit is displayed if the user changes to [°F] in the set-up.

## 6.2 Description of acoustic signals

The acoustic signal provides useful information on key selection and operating processes.

1. A single tone is heard when a key is pressed or the furnace is in "HOLD" mode.
2. A double tone is heard when an invalid key is pressed or invalid data are selected.
3. A five-fold double tone is heard at the end of a combustion cycle.
4. A five-fold single tone is heard if the lift does not reach its final position. Reset by means of the (Lift) key.
5. A continuous uninterrupted tone is heard if the vacuum level cannot be reduced on time. Reset by switching off.

## 7 Putting the unit into operation

### 7.1 Transport

Transport the unit carefully in a horizontal position. Avoid jerking it! For sizes and weight, see Chapitre 11 Technical data

### 7.2 Installation

**CAUTION!**

**Check the contents of the shipping box against the scope of delivery given in the operating manual.**

**See chapter 4 Scope of delivery page GB 38**

Remove all articles from the upper foam section. Carefully remove the upper foam section. Keep the box and packaging for possible later use.

Carefully lift out the furnace, holding it below the muffles on the inside, and place it down upright. Do not lift the furnace by the lift rod.

#### 7.2.1 Vacuum pump connection

1. Before initial operation, the pump must be filled with the special oil provided. Before each further use, the oil level must be checked and topped up if necessary. (See pump operating instructions)
2. Connect the vacuum pump grommet to the rubber vacuum hose.
3. Connect the furnace grommet to the other end of the rubber hose. Make sure that the clamps are secure.
4. Connect the vacuum pump power cord to the mains connection socket on the rear of the furnace. Check the position of the pump: the vacuum hose must not be kinked.
5. Make sure the pump is switched on. For further information, refer to the pump operating instructions.

#### 7.2.2 Furnace connection to the mains power supply

Make sure that the master switch on the rear of the furnace is switched off. Connect the power cord to the furnace mains power connection. Insert the power cord into a wall socket. Switch on at the master switch. A single tone can be heard. Press the "ON/OFF" key at the front (operating panel). The furnace is now in idling mode. Press the [LIFT] key - the platform is lowered automatically. **Place the ceramic table on the platform.** The furnace is now ready for initial operation.

### 7.3 Initial operation


- Before combustion processes are carried out, the furnace must first be subjected to an initial drying process in order to eliminate any moisture that may have penetrated the furnace during transportation.
- The furnace combustion process is a drying process for the combustion chamber and must be carried out before the furnace is put into use.
- The process is described before the operating instructions in order to remind you to read the instructions completely and only then to carry out the initial combustion process before you work with the programs.

#### 7.3.1 Initial combustion process of the furnace

1. Connect the furnace.
2. Lower the lift platform to the lowest position using the "LIFT" key.
3. Place the combustion table on the lift platform.
4. Enter the following parameters in program 99:

START TEMP	500°C
DRY	3 Minutes
Heating speed	40°C / Minute
HIGH TEMP	900°C
HOLD (holding time)	3 Minutes
Without vacuum	

Start program 99. This program eliminates all moisture absorbed in the muffle. Repeat the above-described procedure three to four times.

**NOTE!:**

***The furnace is ready for operation only after this procedure is completed!***

### 7.4 Operating the machine

#### 7.4.1 Switching on

In order to switch on the furnace, move the master switch at the rear to position "I"; the furnace is now in idling mode.

#### 7.4.2 Program mode

The Heramat C ceramic furnace has 100 programs. Programs 0-2 are fixed programs and cannot be altered. Programs 3-29 contain specified data, which can, however, be altered. In programs 4,16,18,29 and 95-98, a dwell period of 4:00 minutes is specified in the lower lift position. This period cannot be altered or deleted. For entering parameter data, see chapter 7.4.7 Entering and altering parameter data page GB 46.

#### 7.4.3 Set-up option


To enter the "Set-up" mode, switch on the furnace (switch "I" on rear) and then select code "50". To select between the options, press the "STEP" key.

1. Press key "1" to select between °C [CELS] and °F [FAHR].
2. Press key "2" to select between mmHg [HG] and mBar [BAR].

3. Press key "3" to select the possibility of programming the vacuum level or use the preprogrammed vacuum level 960mbar/720mmHG (vac. level) as standard: >nvac<: the furnace uses the standard vacuum; >dvac<: activates the manual vacuum programming.
4. Press key "4" to select between ES=Energy Safer active and no ES = Energy Safer deactivated.  
With ES (Energy Safer active), the ceramic furnace automatically switches to stand-by mode after one hour. "ES" appears in the display as program. The combustion chamber closes and the furnace maintains a constant temperature of 200°C. To terminate the Energy Safer mode, the "LIFT UP/DOWN" key must be actuated.
5. Press key "5" to reload original programs. The customer programs 3-29 are deleted, the original programs stored in the E-Prom are read in. When key "5" is pressed, "rElo" flashes in the display and an acoustic signal is given.
6. The software version can now be read out in the last menu item (without number). To leave the "Set-up" mode, press the "on/off" key.

#### 7.4.4 Idling mode

Press the "ON/OFF" key at the front (operating panel); the furnace is now ready for program selection.


**NOTE!:**

*The idling mode is intended for use within extended periods in which the furnace is not in operation. This mode has no influence on the combustion chamber, in which a constant temperature of 300°C is maintained. The lift platform remains in its previous position. It is especially important to keep the combustion chamber closed at all times to prevent moisture from penetrating the muffles. To ensure this, the possibility of switching off the furnace while the lift platform is at the lowest position has been deactivated.*

#### 7.4.5 Program selection

- Press the "PROG" key.
- Use the number keys to enter a number from 30-99 for customer-specific programs or 3-29 for original programs.
- Press the "ENTER" key.

**NOTE!:**

**Please note the following:**

- If you want to change to another program, press the "CLEAR" key and enter a different program number. Confirm with "ENTER".
- During the program mode, the combustion chamber temperature is displayed.
- When the muffle has cooled (below 100°C), >COLD< appears in the display.
- The combustion chamber temperature is maintained at the starting temperature of the selected program (START TEMP).

### 7.4.6 Programming mode

The programming mode serves for entering, checking and altering the data of the program parameters. When you press the "STEP" key, you automatically go to the parameter field, where data can be entered. Press the "PROG" key to leave the parameter field.

### 7.4.7 Entering and altering parameter data

Press the "STEP" keys to go to the desired parameter, and you can read off the value in the data display.

Use the data keys to enter the necessary parameter data. With this procedure, you can move through your parameters and check and alter data as required.


Press the "ENTER" or "STEP" key when you have completed entering or altering the required parameters.

#### Example 1: Entering a START TEMP of 650°C:

1. Press the "STEP" key until the "START" LED lights up.
2. Press the keys "6", "5", "0".
3. Press the "ENTER" or "STEP" key again if you wish to enter other parameter data.

#### Example 2: Entering a VACUUM HOLD period of 1 min. 35 s.

1. Press the arrow key until the "HOLD" LED lights up.
2. Press the keys "1", "3", "5".
3. Press the "ENTER" or "STEP" key again.

**NOTE!:**

*The (Minutes : Seconds) display stops flashing after confirmation. (If a Vacuum Hold period is entered, the V Off value disappears).*

### 7.4.8 Description of parameters and their limit ranges

Parameter	Function	Range	Remarks
START	The start temperature of the program. The combustion process begins as soon as it is reached.	300°C to 800°C	One of the parameters that must be entered to run the program.
DRY	The time in which the lift platform moves out of the lowest position until the muffles are closed.	0 to 30 minutes	While DRY is displayed, the lift platform moves up in five steps. The original programs 4, 16, 18+29 (paste-opaque) use a waiting period of 4 minutes in lowest lift position, in addition to the Dry time. The customer programs 95-98 also use this waiting period of 4 minutes in addition to the freely programmable Dry time.
PRE HEAT	The final stage of DRY, during which the muffle is closed and the temperature does not rise for the specified period.	0 to 30 minutes	The vacuum pump starts up only after the PRE-HEAT period has expired.

Parameter	Function	Range	Remarks
HEAT RATE	The speed in °C or °F per minute at which the temperature rises from START TEMP to HIGH TEMP and, if applicable, from HIGH TEMP to TEMPER (hardening).	10°C to 180°C (17°F to 324°F)	One of the parameters that must be entered to run the program.
HIGH TEMP	The temperature to be reached at the end of HEAT RATE, or the intermediate temperature if the value of TEMPER is greater.	300°C to 1200°C (572°F to 2,192°F)	Must not be below "VACUUM OFF" temperature. (If a vacuum is specified in the program). One of the parameters that must be entered to run the program.
HOLD	The period in which the furnace remains at HIGH TEMP before the temperature rises/falls to TEMPER (if specified).	0 to 30 minutes (stated in minutes and seconds)	Intermediate delay if TEMPER is greater than HIGH TEMP.
TEMPER	An additional temperature stage. Can be lower or higher than HIGH TEMP. If higher, heating process is restarted. If lower, chamber is opened a little to allow rapid cooling to TEMPER value. The lift platform then closes the muffle and the temperature is stabilised.	100°C - 1200°C (212°F to 2,192°F)	
HOLD	The period in which the furnace remains at TEMPER before the chamber is opened.	0 to 30 minutes	Only use if TEMPER parameter has been programmed; if necessary, delete values with "CLEAR" [ _ _ ], do not press "0".
COOL TIME	The time required to lower the lift platform. The lift platform is lowered in four consecutive steps, starting with the closed combustion chamber to the lowest position. Heating is switched off.	0 to 30 minutes	
V. ON (Vacuum on)	Vacuum starts with the specified temperature.	Each value within the selected temperature range of the program	If the value is the same as START TEMP, the vacuum is created after PRE-HEAT has expired.
V. OFF (Vacuum off)	The temperature at which the vacuum is ended.	Each value within the selected temperature range of the program	Must not be greater than HIGH TEMP or TEMPER (if TEMPER has been set higher than HIGH TEMP). Only (VAC.) OFF or (VAC.) HOLD may be entered.

Parameter	Function	Range	Remarks
HOLD (Hold vacuum)	Period for maintaining the vacuum with HIGH TEMP or TEMPER (if TEMP HOLD has been specified)	0 to 30 minutes	(VAC.) HOLD must not be greater than TEMP HOLD. Only (VAC.) OFF or (VAC.) HOLD may be entered.
VAC (Vacuum level)	The level at which the vacuum is maintained during the combustion process.	From 100mbar to 980 mbar or 100 mmHg to 760 mmHg >dvac< at 990 mbar continuous operation >dvac< at 760 mmHg continuous operation	In the standard setting, the value is specified and not editable (960 mbar / 720 mmHG). For editing, the setting "dVAC" must be activated in the setup (see chapter 7.4.3 Set-up option page GB 44). For a combustion process with vacuum, VAC.ON and VAC.OFF or HOLD must always also be programmed.

**NOTE!:**

*If the TEMPER parameter has been specified, and the value is greater than HIGH.TEMP, the vacuum is maintained with TEMPER for a specified (VAC.) HOLD (hold vacuum) period. If there is no TEMPER value or TEMPER is lower than the HIGH TEMP value, (VAC.) HOLD (hold vacuum) depends on HIGH TEMP HOLD (hold temperature).*



Fig. 6 Diagram


**NOTE!:**  
**Please note the following:**

1. If you program these data in a parameter, check that the data are correct and the correlation with the other parameters.
2. Please take into consideration that parameters have limit ranges; if you enter data above or below these limits, the data will not be accepted (double tone).
3. In the event of a fault, the warning tone sounds when the "START/STOP" key is pressed, and the combustion process is not initiated.

## 8 The combustion process

### 8.1 Starting the combustion process

1. Make sure that the lift platform is at the lowest position. Briefly press the LIFT key to lower the lift platform, press again to raise the lift platform. The combustion process will not be started until the lift platform has been completely lowered.
2. Select a program by means of the procedure described. (See chapter 7.4.2 Program mode page GB 44.)
3. Press the "START/STOP" key.
4. The combustion process starts automatically as soon as the combustion chamber temperature reaches the temperature value for START TEMP and the relevant function LED lights up.
5. When the combustion process is completed and the lift platform has been lowered to the lowest position, a five-fold double tone is emitted and the red END lamp lights up.


**NOTE!:**  
**Please note the following:**

If the program includes a vacuum, the vacuum lamp lights up at the start of pump operation and the microprocessor monitors the function during the entire procedure.

- If the vacuum pump is not connected or is not running, VAC lights up and a double warning tone is emitted. The combustion process is interrupted and the platform is lowered. VAC starts flashing. Eliminate the fault and restart the program.
- Only when the set vacuum level is reached does the temperature start to rise (this is particularly useful for ceramics which are burnt for a brief period).
- Dental work must only be placed on the tray. Avoid placing objects or other work on other areas of the furnace.

## 8.2 Idling mode

We recommend that the furnace be kept in idling mode with the muffle closed until you wish to start the next combustion process.

**NOTE!:**

*It is not possible to switch the furnace off (On / OFF) if the lift platform is not in the uppermost (closed) position. In order to switch off the furnace, please close the chamber using the LIFT key.*

## 8.3 Interrupting the combustion process

If you wish to interrupt a program during the combustion process, please press the "START/STOP" key. The vacuum in the chamber is reduced, the lift platform is lowered and the furnace remains in the idling mode.

## 9 Additional features

### 9.1 Time remaining to end of combustion process

You can read off the remaining time by pressing key "1".

### 9.2 Post Fire

Towards the end of the combustion process during the cooling period, before the END display appears, you have the chance of post-firing the work at a higher temperature. You can use the "HOLD" key (9) to hold the lift platform (27) in the lower position.

Press the "P. FIRE" key in the post-firing process to move the lift platform back into the combustion chamber and raise the temperature by 5°C.

Each subsequent immediate pressing of the "P. FIRE" key raises the temperature by an additional 3°C.

All other parameters, e.g. holding time, remain unaltered. At the end of the combustion process, the maximum temperature is automatically reset to the original programmed value.

### 9.3 The HOLD key


When the "HOLD" key is pressed, the combustion process is temporarily interrupted. The temperature is maintained and the lift platform moves no further.

### 9.4 Altering parameters during the combustion process

In order to alter parameters during the combustion process:

1. Press the "HOLD" key.
2. Go to the parameter you wish to alter (using the "STEP" key).
3. Make the required alteration.
4. Press the "HOLD" key again.

The furnace carries out the combustion process according to the new parameters. The parameters are not stored, but reset after completion of the program to the original value (especially important for the post-firing function!).

**NOTE!:**  
**Please note the following:**

- Press once to put the furnace in a "HOLD" phase.
- Press again to make the furnace end the "HOLD" phase (the combustion process is continued).
- During the "HOLD" phase, a single tone sounds every ten seconds to remind you that the furnace is in the "HOLD" phase.
- If you press the "HOLD" key after pressing the "LIFT" key, the LIFT movement is stopped.

## 9.5 Monitoring the vacuum level

If key "2" is pressed during the combustion process, it is possible to monitor the vacuum level on the data display. After 5 seconds, the display automatically reverts to the temperature.

## 9.6 Night mode

There are two ways to enter the night mode:

Either press the "ON/OFF" key during a combustion process. After completing the combustion process, the furnace automatically enters the night mode. The automatic night mode is indicated by the flashing of the combustion process indicator lamp located in the "START/STOP" key.

The night mode comes into operation only after the combustion chamber temperature has dropped to 100°C; it then causes the platform to close and maintains the chamber at a constant temperature of 100°C.

In order to interrupt the automatic night mode, again press the "ON/OFF" key (1).

You can start the manual night program "0" by pressing "PROG", "0", "ENTER", "LIFT". The lift closes the muffle and the temperature drops during the night mode to 100°C. To interrupt the night mode, press the (STOP) key.

## 9.7 LIFT key

The "LIFT" key serves to move the lift platform. Briefly press the "LIFT" key to lower the lift platform, press again to extend the lift platform.

## 9.8 STOP key

If you wish to interrupt a program during a combustion process, please press the "START/STOP" key. The vacuum in the chamber is reduced, the lift platform is lowered and the furnace remains in the idling mode.

## 9.9 Temper

If the TEMPER (harden) parameter is lower than HIGH TEMP, TEMPER and TEMP HOLD act as a cooling stage, i.e. the temperature drops, with the muffle opened slightly, to the specified TEMPER temperature and the muffle closes again. The temperature remains at the value for the specified TEMP HOLD period.

These two special parameters, which have been integrated in the HERAMAT C ceramic furnace, serve to adjust the combustion process to the expansion coefficients, both for metal and ceramic, and for the first time provide flexibility in heating processes for future types of ceramic.

## 9.10 Diagnosis mode

To enter the diagnosis mode, the furnace must be switched on (switch "1" on rear); then enter the code "562". This mode allows the checking of various functions of the furnace, e.g.:

1. Checking the muffle heating (by pressing key "1"). You can see how the temperature rises by pressing key "1".
2. Checking the vacuum system (by pressing key "2"). The pump starts to work and the vacuum level rises). When you release the key, the vacuum remains at the level reached. If you press again, the vacuum is reduced.
3. Checking the Lift UP movement (by pressing key "3").
4. Checking the Lift DOWN movement (by pressing key "4").
5. If you enter the code "562042", you can see the entire operating hours of the furnace above 600°C.

To select between these options, press the "STEP" key (up - down).

To leave the diagnosis menu, actuate the mains switch "I/O".


## 10 Servicing

### **CAUTION!**

*The functional efficiency and safety of the unit is only guaranteed if the necessary tests, service, and repair work is performed by trained service personnel. Heraeus Kulzer GmbH will accept no liability if the unit is damaged by faulty repair work which was not performed by trained service personnel, or if original spare parts and accessories are not used to replace parts.*

You are advised to conclude a service contract with our service centres - you can request an offer from them (addresses: see chapter 13 Heraeus Kulzer addresses worldwide page GB 60).


### **ELECTRIC SHOCK!**

*Before using the unit, check the mains lead and plug for damage. If they are damaged, do not connect the unit to the mains.*


Wipe down the outside and control elements in the unit with a mild soap solution (rinsing solution) and a soft, damp (but not wet!) cloth. Check that the vacuum pump functions perfectly with the combustion chamber empty every 2 weeks. Watch the display (key "2"). The display should have reached the set final value (e.g. 960 mbar) within max. 30 seconds after the vacuum pump is switched on.

If the display does not indicate this, or is significantly slower:

- The gaskets on the lift or combustion chamber may be dirty or damaged
- The vacuum pump may not be working properly (check oil)
- The limit switch may be incorrectly set (lift does not close completely). Poor combustion may result. Please contact your Service partner (see chapter 13 Heraeus Kulzer addresses worldwide page GB 60).


**NOTE!:**  
Please note the following:

**DISCONNECT FROM MAINS!** Before carrying out maintenance work on the furnace, or in event of a mechanical failure, switch off the power supply (mains switch "I/O") at the rear of the furnace and disconnect the unit from the power supply (disconnect mains connector).

- The furnace must be opened only by a service technician authorised by the manufacturer or the furnace supplier. There are no components requiring maintenance by the customer inside the furnace.
- Unauthorised opening of the furnace will lead to cancellation of warranty.
- Pumps connected to the furnace by the manufacturer or another supplier must be designed and suitable for the intended use of the furnace, and bear the CE mark. In event of any doubt, please contact Heraeus Kulzer GmbH, Hanau.
- If it is necessary to change the fuse:

**DISCONNECT FROM MAINS!** Danger of electric shock if unit is opened  
Unplug the unit before opening it.

1. Disconnect the unit from the mains supply.
  2. Replace the fuse with the spare fuse provided or a fuses of equivalent quality.
  3. Check whether the spare fuse is of the correct value. Fuses with deviating data must not be used. Danger to property and persons!
  4. In event of repeated failure, take the unit out of operation (pull mains connector, secure against operation) and inform the relevant Service.
- In event of an extended external power failure while the furnace is in operation: disconnect the unit from the mains power supply. Lower the table by inserting the wrench provided into the hexagon screw at the rear of the furnace.
  - Further working is not permissible and must therefore not be carried out.

## 10.1 Tests

Check the following features in the unit at least once annually, and ensure it is in working condition:

- Mechanics
- Function (technical data)
- Electrics
- Safety devices and features of the unit (e.g. signal tone, lift movement, start/stop)
- Servicing

Work on the unit's electrical equipment must only be performed by trained service personnel. The unit must be switched off when such work is performed. Only use authorised original spare parts.

## 10.2 Preheating

We recommend that the furnace be preheated to a working temperature before you carry out calibration measures. We suggest a preheating procedure in which you carry out three times consecutively a combustion process with parameters similar to those in the following Table 10. After completion of the combustion processes, leave the furnace for 30 minutes at 500°C.

- Move the lift to the lower position. (LIFT key).
- Connect the vacuum gauge between the furnace and the vacuum pump. (Use T connection).
- Insert the connector into the calibration socket at the rear of the furnace.

## 10.3 Vacuum calibration

In order to start the calibration process, press consecutively the following key: "PROG", "No.1", "ENTER" and "START".

The lift platform rises step by step to the closed position.

When the lift platform is in the fully closed position, press key "8" to switch on the pump.

Observe the vacuum gauge. AS soon as the gauge reaches 960 mbar or 720 mmHg, press key "1"; a double tone is emitted and the pump stops. Press Stop. The lift platform is lowered.

Only after the table has been moved to the lowest position may the calibration connector be pulled from the socket at the rear and the vacuum gauge separated from the vacuum hose. Reconnect the pump correctly.

The furnace is now ready for operation.

**Tabel 1. Heramat C preset programs - New CPU with 100 (97) variable programs - from software V 2.00 on**

	Heramat C	STAR T	DRY	PRE HEAT	HEAT RATE	HIGH TEMP	HOLD	TEMPER	TEMP HOLD	COOL TIME	V ON	V OFF	V HOLD
		[°C]	[min]	[min]	[°C/min]	[°C/min]	[min]	[°C]	[min]	[min]	[°C]	[°C]	[min]
<b>Furnace function</b>													
0	Night mode - not alterable	100											
1	Silver-test - not alterable	500	-	-	50	1060	-	-	-	-	1059	1060	-
2	Furnace cleaning - not alterable	600	1:00	1:00	100	1100	20:00	-	-	-	600	-	20:00
<b>HeraCeram</b>													
3	Oxide 880	600	-	-	100	880	10:00	-	-	-	600	-	10:00
4	Paste-opaque special program*	600	*2:00	1:00	100	880	1:00	-	-	-	600	880	-
5	Paste-opaque normal program	600	5:00	1:00	100	880	1:00	-	-	-	600	880	-
6	Powder opaquer	600	2:00	1:00	100	880	1:00	-	-	-	600	880	-
7	Silver-test high-fusing margin	600	3:00	1:00	100	870	1:00	-	-	-	600	870	-
8	Shoulder-material high-fusing margin	600	2:00	1:00	100	860	1:00	-	-	-	600	860	-
9	dentin basic firing cycle	600	3:00	2:00	100	860	1:00	-	-	-	600	860	-
10	Dentin corrections	600	3:00	2:00	100	850	1:00	-	-	-	600	850	-
11	Gloss 850	600	2:00	2:00	100	850	0:30	-	-	-	-	-	-
12	Gloss 840	600	2:00	2:00	100	840	0:30	-	-	-	-	-	-
13	Correction compound	600	2:00	2:00	100	810	1:00	-	-	-	600	810	-
14	Shoulder-material low-fusing margin	600	2:00	2:00	100	790	1:00	-	-	-	600	790	-
15	Soldering (recommended solders - see below)	600	5:00	5:00	55	850	-	-	-	-	-	-	-
16	Blendgold new HeraCeram Preciano Special pr.* Blendgold HeraCeramSUN	400	*2:00	-	55	900	-	-	-	-	-	-	-
Special program* = 4:00 min. in lowest lift position + 2:00 min. furnace closed time = Specified drying time = furnace closed time													
Recommended solders: Aura Solder 800, Herador Solder 800, Herador PF Solder 800, Maingold Solder S 800, Maingold Solder 800, Maingold PF Solder 800													

**Tabel 2. Heramat C preset programs - New CPU with 100 (97) variable programs - from software V 2.00 on**

Heramat C		STAR T	DRY	PRE HEAT	HEAT RATE	HIGH TEMP	HOLD	TEMPER	TEMP HOLD	COOL TIME	V ON	V OFF	V HOLD
		[°C]	[min]	[min]	[°C/min]	[°C/min]	[min]	[°C]	[min]	[min]	[°C]	[°C]	[min]
<b>HeraCeram Sun</b>													
17	HeraSun Oxide	600	1:00	-	100	800	10:00	-	-	-	-	-	-
18	Basic / paste-opaque Special Progr.*	600	*2:00	1:00	100	790	1:00	-	-	-	600	790	-
19	Basic / paste-opaque Normal Progr.	600	5:00	1:00	100	790	1:00	-	-	-	600	790	-
20	Powder opaquer	600	3:00	1:00	100	790	1:00	-	-	-	600	790	-
21	Shoulder-material high-fusing margin	600	3:00	1:00	100	775	1:00	-	-	-	600	775	-
22	dentin basic firing cycle	600	3:00	2:00	100	760	1:00	-	-	-	600	760	-
23	Dentin corrections	600	3:00	2:00	100	750	1:00	-	-	-	600	750	-
24	Gloss 750	600	3:00	1:00	100	750	0:30	-	-	-	-	-	-
25	Gloss 740	600	3:00	1:00	100	740	0:30	-	-	-	-	-	-
26	Correction compound	600	3:00	1:00	100	710	1:00	-	-	-	600	710	-
27	Shoulder-material low-fusing margin	600	3:00	1:00	100	700	1:00	-	-	-	600	700	-
28	Soldering HeraCeram Sun Solder 2	600	8:00	5:00	55	740	0:30	-	-	-	-	-	-
29	Blendgold New Special Program*	400	*2:00	-	55	820	-	-	-	-	-	-	-
Special program* = 4:00 min. in lowest lift position + 2:00 min. furnace closed time = Specified drying time = furnace closed time													
Recommended solders: Aura Solder 800, Herador Solder 800, Herador PF Solder 800, Maingold Solder S 800, Maingold Solder 800, Maingold PF Solder 800													

## 11 Technical data

### 11.1 General

Preset combustion programs	30
User-defined combustion programs	70
Dimensions	B 30 cm x T 33 cm x H 58 cm W 12" x D 13" x H 23"
Muffle chamber	dia. 95 mm x 60 mm
Combustion table	dia. 90 mm
Net weight	22 kg ( 48.5 lbs )
Mains voltage	230 V 50/60 Hz – 1800 W 115 V 50/60 Hz – 1400 W
Pump capacity	230 V / 5 AMP
Temperature range	100°C - 1200°C
Heating speed Heating speed	230 V from 10°C/min. to 180°C/min 115V from 10°C/min. to 120°C/min
Temperature preakup	1°C
Vacuum level	0 to -990 mbar under atmospheric pressure or 100 to 750 mm Hg under atmospheric pressure
Raising/lowering time for table	20 seconds - 30 minutes
Temperature delay	0 seconds - 30 minutes
Vacuum delay	0 seconds - 30 minutes
Time accuracy (steps)	1 second

### 11.2 Transportation and storing conditions:

Temperature 0°C to + 40 °C  
Relative air humidity less than 100% rel. moisture

### 11.3 Operating conditions:

Temperature 10°C to + 30 °C  
Air humidity 80% rel. humidity up to 31°C  
50% rel. humidity at 40°C

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**Fig. 7 Heramat C rating plate**


## 11.5 Disposal of old equipment according to WEEE

### Electrical and Electronic Equipment Act (ElektroG)

This Act sets out requirements for electrical and electronic equipment pursuant to directive 2002/96/EG issued by the European Parliament and the European Council of 2005-05-03. Its main purpose is to prevent waste from electrical and electronic equipment and to promote reuse, recycling and other forms of recovery to reduce both the volume of waste for disposal and the inclusion in waste of harmful substances from electrical and electronic equipment.

Mandatory labelling will be effective from August 13th, 2005 and includes products which were produced and delivered after this date.

**Consequently, corresponding products from Heraeus Kulzer will be marked with the symbol below.**

For detailed information regarding the proper disposal of equipment which is no longer in use, please read the operating instructions of the unit or contact your local dealer.


### **IMPORTANT!**

**Marked equipment must not be brought to local waste disposal centers.**

## 12 Circuit diagrams

### 12.1 Terminal allocation scheme

## 13 Heraeus Kulzer adresses worldwide

COUNTRY	NAME ADRESS
<b>Deutschland/Germany</b>	<b>Heraeus Kulzer GmbH</b> Produktsupport Labor Grüner Weg 11 63450 Hanau Tel. +49 (0) 6181 / 35-5894 oder 35-4773 Fax +49 (0) 6181 / 35-5993
<b>Australien/Australia</b>	<b>Heraeus Kulzer Australia Pty. Ltd.</b> Unit 18A / 390 / Eastern Valley Way Roseville NSW 2069 Locked Bag 750 Roseville NSW 2069 Tel. +61 294178411 Fax +61 294175093
<b>Brasilien/Brazil (America Sul/ America del Sur)</b>	<b>Heraeus Kulzer South America Ltda.</b> Rua Lisboa, 90 005413-000 Sao Paulo SP Tel. +55 11 30688171 Fax +55 1130688172
<b>China/China</b>	<b>Heraeus Kulzer Dental Ltd.</b> 680 Guiping Road 200233 Shanghai Tel. +86 2164958488 Fax +86 2164951732
<b>Frankreich/France</b>	<b>Heraeus S.A.S.</b> Parc Silic - Bat.i.2 Villebon - BP 630 12, Avenue du Québec 91945 Courtaboeuf Cédex Tel. +33 169184848 Fax +33 169287822
<b>Großbritannien/United Kingdom</b>	<b>Heraeus Kulzer Ltd.</b> Heraeus House Albert Road / Northbrook Street RG14 1DL Newbury, Berkshire Tel. +44 1 63530500 Fax +44 1 635524622
<b>Indien/India</b>	<b>Heraeus Kulzer Dental India Private Ltd.</b> Heraeus Kulzer House 344/2 Ladoo Sarai 110030 New Delhi Tel. +91 116512849 Fax +91 116512869
<b>Italien/Italy</b>	<b>Heraeus Kulzer S.r.l.</b> Via Console Flaminio 5/7 20134 Milano Tel. +39 022100941 Fax +39 02 210094283

COUNTRY	NAME ADDRESS
Japan/Japan	<b>Heraeus Kulzer Japan Co. Ltd.</b> 9F Crystal Tower 2-27 Shiromi 540-6009 Osaka Tel. +81 64791 3550 Fax +81 647913561
Mexiko/Mexico	<b>Heraeus Kulzer Mexico S.A. de C.V.</b> Homero 527 - 301 y 302, Col. Pol. 11560 Mexico Tel. +52 5 5315549 Fax +52 5 552551651
Niederlande/Netherlands	<b>CAVEX Holland B.V.</b> Harmenjansweg 19 - 22 2003 RW Haarlem Tel. +31 23-5307700 Fax +31 23-5356482
Nordamerika/North America	<b>Heraeus Kulzer Inc.</b> Headquarter 99 Business Park Drive 10504 Armonk Tel. +1 914 273 8600 Fax +1 914 273 5211
Österreich/Austria	<b>Heraeus Kulzer Austria GmbH</b> Kaiserstraße 119 1070 Wien Tel. +43 (0) 1 4080941 Fax +43 (0) 1 4080941-70
Skandinavien/Scandinavia	<b>Heraeus Kulzer Nordic AB</b> Hammarbacken 4 19149 SOLLENTUNA Tel. +46 85 8577755 Fax +46 86 231413
Spanien/Spain	<b>Heraeus S.A.</b> Manuel Tovar 24 28034 Madrid Tel. +34 913580375 Fax +34 913580368

## 14 Document history

- 26.03.03 First edition
- 07.06.05 Company name and name plate updated
- 04.09.06 Added chapter 11.5 Disposal of old equipment according to WEEE

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We reserve the right to make technical changes!

**Heraeus Kulzer GmbH**

Grüner Weg 11  
63450 Hanau

Telefon: +49 6181 35-5894

Telefax: +49 6181 35-5993

info@heraeus-kulzer.com

www.heraeus-kulzer.com