
## Palajet®

- Verfahrens- und Gebrauchsanweisung
- Method and user instruction
- Methodé et mode d'emploi
- Procedimento ed istruzioni di impiego
- Modo de empleo del proceso y uncionamiento

# Palajet®

Injektionsverfahren zur Herstellung von Prothesen  
Injection process for the fabrication of dentures  
Système d'injection pour la fabrication de prothèses  
Procedimento di iniezione per la realizzazione di protesi  
Método de inyección para la fabricación de prótesis

## Bedienungselemente

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## Palajet®

### Area of Operation

Palajet®, air pressure injection apparatus for the production of dentures with cold polymer PalaXpress® or hot polymers such as Paladon 65® and Palaimpact®.

### Equipment Characteristics and Process

The Palajet is a pneumatic injection apparatus for injecting denture material. The air inlet should not exceed a maximum of 10 bar. This apparatus is set to a working pressure of 4 bar, which can be controlled with the pressure meter. The plunger speed during injection is 8-10 mm/s. The Palajet® is designed and developed according to the most common safety regulations and guidelines.

The injection procedure runs automatically after the switch has been flipped. The concluding polymerisation of the dentures occurs in the polymerisation apparatus Palamat® elite.

During the injection process the mixed acrylic dough is pressed through the injection cylinder into the denture flask cavity with the help of the pressure plunger. After the filling process is completed, the remaining denture acrylic will escape through the ventilation channel which is then closed. Please observe the deviation when using hot-melt polymers! The complete polymerisation of the denture occurs in the polymerisation apparatus Palamat® elite.

■ EC – Declaration of conformity  
We, Heraeus Kulzer GmbH, Grüner Weg 11, 63450 Hanau, hereby declare that the machine named below is designed and built to meet the basic relevant health and safety requirements specified in the EC Directive, as is the model launched by us.

If the machine is changed without our approval, this declaration will cease to be valid.

Relevant EC directive: Machine Directive 98/37 EG, Appendix I

### Starting Apparatus

This apparatus must be positioned in a dry room only and on a solid foundation so that it cannot overturn. Prior to the compressed-air connexion control that the primary pressure is not more than 10 bar and not less than 4 bar. At the same time it must be insured that the air pressure is free of oil and water. The apparatus is connected to the air circuit with a hose (inner Ø 8 mm).

#### ■ Preparation

After compressed-air connexion flip the switch down so that the pressure plunger releases the fill cylinder (pressure plunger moves down). Move the pressure module to the rear as far as it will go. Remove the fill cylinder from the apparatus by a left-hand rotation and take it out. Unscrew the fill cylinder cover from the fill cylinder.

#### ■ Investment

Before the embedding make sure that both flask halves fit exactly (check for the

same identification number). Lubricate the inside with Vaseline so that the plaster can later be easily removed. The spacer (placeholder) can be used to save plaster if the Duoflask is used with only one model.

The fabricated wax dentures are carefully positioned so that when embedded the injection channel can be attached on the upper jaw at the tuberculum maxillae or the lower jaw at the trigonum retromolar.

Remove excess plaster and smooth the surface. Clean the edges of the flask from excess plaster. Use class III plaster for embedding. After the plaster has hardened, attach the injection channel (Ø 7 mm) and the wax for the ventilation channel (not required for the Paladon® 65 and Palaimpact®).

#### PalaXpress®:

Position the ventilation channel (Ø 3 mm) in the area of the tuberculum maxillae or trigonum retromolar (when using the single flask position the ventilation channel on the opposite side of the injection channel), choosing the shortest distance to the ventilation channel.

Insulate the plaster surface with Aislar®. Place the flask halves on top of each other and firmly screw them together (without using a hammer). Through the flask opening fill in plaster so that the teeth are covered and let harden. Insulate the plaster surface with Aislar® and fill flask with plaster. Wipe off remaining material so that the flask fits between the centring disks of the Palajet®.

#### ■ Insulating

After the plaster has hardened, preheat the flask in hot water (approximately 90° C) without chemical additives in order to soften the wax, separate the flask halves after circa 5 minutes and remove the remaining wax with hot, clear water. Thinly insulate the hot plaster on the contact area to the denture plastic twice (allowing each layer to dry) with Aislar® (avoid building pools).

Roughen the basal area of the teeth after insulation of the plaster with a coarse diamond (remove sanding dust). To improve the connection between the teeth and the denture plastic, add the adhesive agent Palabond® twice with a brush. Allow each layer to dry for 30 seconds. After the second layer the adhesive agent will remain active for 10 minutes.

Close both flask halves firmly and insert the flask into the Palajet®. The sealing surface of the vent seal of the Palajet apparatus should not abut against the opening of the flask's ventilation channel.

#### ■ Preparing acrylic

Prior to the mixing procedure unscrew the cylinder cover from the fill cylinder and insert the fill cylinder cartridge. In order to receive the acrylic dough, the seal disk with the sealing lip showing upwards must be pressed into the fill cylinder cartridge (as far as it will go). The cover cartridge will be inserted into the cylinder cover.

Pour the powder continuously into the reserved liquid, mix well for 30 seconds (transfer acrylic dough to the filling cylinder immediately) and cover.

	Mixing ratio
PalaXpress®	2 g : 1 ml 30 g : 15 ml
Paladon®65	5 g : 2 ml 30 g : 12 ml
Palaimpact®	21 g : 10 ml 32 g : 15 ml

### Note

Do not wipe off the mixing cup on the surface of the fill cylinder cartridge to avoid inhomogeneities.

Mix in one rotation to avoid incorporating air bubbles. When working with PalaXpress® pour the acrylic dough in a thin jet free of bubbles into the prepared fill cylinder cartridge.

#### ■ Injecting

PalaXpress®:

After filling the fill cylinder cartridge wait for the moment of injection. This moment is reached when the dough maintains a dull surface. Please observe the deviation when using hot-melt polymers.

Paladon® 65 / Palaimpact®:

From the start of the mixing process, the precise waiting period before injection, at 23°C, is as follows:

Paladon® 65 – 20 minutes

Palaimpact® – 10 minutes.

Waiting times are dependent on the room temperature and the amount that was mixed.

Close the fill cylinder with the cylinder cover and turn it into the Palajet®. Lock the flask in position with the flask holder knob.

Paladon/Palaimpact:

Close the air drainage channel with the device closure screw prior to injection.

After pulling the pressure module as far as it will go, the switch is reversed upwards while the fill cylinder is inserted and the denture plastic is injected with help of the pressure plunger. As soon as material escapes from the opening of the ventilation channel, the ventilation seal must be closed.

PalaXpress®:

Leave the flask 5 minutes in the Palajet® until the acrylic dough no longer absorbs pressure, which means until no more acrylic dough can escape when opening the ventilation seal. To control the time, set timer to 5 minutes.

Paladon® 65 / Palaimpact®:

Allow the flask to remain in the Palajet®, for 25 minutes

#### ■ Polymerising

For the polymerisation of the material the flask must be removed from the injection apparatus. Place the switch in the starting position; the pressure plunger exits the injection cylinder. Push the pressure module to the rear as far as it will go. After screwing out the fill cylinder the flask can be removed, placed on the flask filling basket and therefore inserted into the polymerisation unit Palamat® elite.

PalaXpress®:

The time required for polymerisation is 30 min at 55°C.

Paladon® 65 / Palaimpact®:

Rapid polymerisation takes 30 min at 90°C.

Polymerisation in stages takes 2 h at 70°C – 90°C.

Slow polymerisation takes 10 h at 90°C

Turn back the vent seal and remove any acrylic remains. Fill cylinder and cover cartridge are constructed for multiple use (approximately 5 times). For better sealing and a secure reception of pressure a new sealing disk must be used for each injection procedure!

The removal of remaining acrylic from the fill cylinder and the cover cartridge is done after hardening is completed.

Take care that the cartridges are not damaged!

## Cleaning and Maintenance

Clean the surface of the apparatus with a moist cloth. Make sure that no water can enter.

Do not use sharp cleaning agents. Control the water separator weekly and eventually remove collected water. A pressure control is recommended every two year.

## Safety Notes

Maximum hoisting speed 10 mm/s. Flip the switch upwards only when the fill cylinder is inserted. Before working on the pressure module remove the apparatus from the pressure supply!

### ■ Deflasking

After the polymerisation, the flask is removed together with the flask filling basket from the Palamat® elite. Prior to deflasking the dentures let the flask cool to room temperature. Afterwards follows the deflasking of the flask.

Remove the plaster with a plastic hammer from the flask halves.

### Note

Do not use a metal hammer! When using metal tools on the flask it will be deformed; this can result in inexact fitting precision. The final occlusion must be controlled prior to removing the denture from the model!

### ■ Preparation for further beginning of operation

Prepare flask before reusing and if necessary lubricate with Vaseline. Clean the sealing surfaces as well as the injection and ventilation channel from plaster and plastic remains in order to avoid raising bite.

## Sources of error

Error sources	Remedy
Incorrect operating pressure	<ul style="list-style-type: none"> <li>- Pull setting knob on the regulator- / filter unit upwards until raster snaps out</li> <li>- Adjust regulation of operating pressure to 4 bar and press setting knob downwards until it snaps in</li> </ul>
Too fast or too slow injection speed	<ul style="list-style-type: none"> <li>- Remove protective stirrup above the throttling check valve</li> <li>- Remove locknut on the throttling check valve</li> <li>- Adjust knurled head screw according to the desired speed and fasten the locknut</li> <li>- Hoisting speed = upper throttling check valve</li> <li>- Sink speed = lower throttling check valve</li> </ul>

## Exclusion of Liability

For other usage differing from the instruction of operation and for unauthorized changes no liability will be accepted from the manufacturer. Repairs are allowed only by the manufacturer or by authorized specialized dealers using original parts.

## Technical data

Pressure connection	max. 10 bar
Nominal pressure	4 bar
Piston speed	8 - 10 mm/s
Pressing force on the fill cylinder	approx. 200 kp
Height Width Depth	ca. 400 mm ca. 290 mm ca. 265 mm
Weight	13,5 kg
Apparatus/Serial No.	516

Technical modifications reserved!

### Note

Store the operating instructions visibly near the apparatus!

## Delivery units

### ■ Apparatus Base unit

Palajet®

Injection cylinder including cover

1 flask

12 plastic cylinder cartridges including

12 cover cartridges and

50 sealing disks

Measurement cup

Mixing cup

Mixing spatula

Channel wax (33 strips Ø 7 mm and

30 strips Ø 3 mm)

Timer

Flask filling basket

### ■ Apparat Supplement unit

Palajet®

Injection cylinder including cover

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Directives CEE applicables: Directive machine 98/37 CEE, annexe I

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## Preparation

- Form the wax-up anatomically.
- Embedding of the denture in the base of flask (with out guide pins) with plaster (class III or IV).
- Smoothing of the plaster surface
- Placing of the injection channel (UJ = tuberculum maxillae/LJ = trigonum retromolar ▶ 1
- Isolation of the plaster surface with Aislar®.
- Place the flask halves on top of each other and firmly attach them together (without using a hammer). Fill with plaster so that the dental arch is covered and allow to harden
- Filling flasks with plaster, smooth evenly and allow to harden.
- Place the flasks for 5 minutes in hot water (circa 90°C).
- Open flasks and remove the remaining wax with hot, clear water (without chemical additives
- Isolate the hot plaster surface by applying Aislar® thinly on the contact surface of the denture plastic (avoiding pools).
- Roughen the basal area of the teeth with a coarse diamond (remove dust par-

- tics), apply Palabond® twice and allow each layer to dry 30 seconds. After the second layer Palabond® will remain active for 10 minutes
- Close both flask halves firmly.
- Inserting the flask into the Palajet®.

## Proportioning and mixing

- The denture acrylic is measured and mixed according to the respective manual of instruction. ▶ 2
- Place the seal disk with the sealing lip showing upwards into the fill cylinder cartridge and push as far as it will go.
- Pour the bubble free liquid into the fill cylinder cartridge without touching the edges and avoid wiping off the mixing cup. ▶ 3
- Close the fill cylinder after the moment of injection is reached. ▶ 4
- Screw the fill cylinder into the apparatus.
- Close the flask and tighten screws.

## Injecting

- Pull the pressure module forward and flip the switch upwards.
- Close the vent valve at the ventilation channel after excess material escapes. ▶ 5

## Polymerising

- Leave the flask under injection for 5 minutes in the apparatus, then remove and polymerise according to the manual of instruction in the Palamat® elite. ▶ 6
- Remove flask with the flask filling basket from the Palamat® elite and let cool to room temperature for 30-60 minutes.

## Final handling

- For deflasking of the denture unscrew and open the flask.
- Remove the plaster with a plastic hammer from the flask.
- Examine the final occlusion before removing the denture from the model.
- Remove the denture from the model and begin final finishing.

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## **Heraeus KulzerSouth America Ltd.**
## **Heraeus S.A.**
## **Heraeus Kulzer Dental India PVT Ltd**
### **Pty Ltd**