
Siloc

Verfahrens- und Gebrauchsanweisung
Method and User instruction
Mode d'emploi
Istruzioni per il procedimento e l'impiego
Instructivo de uso y manejo
Instruções de uso

Procedure- en Gebruiksaanwijzing
Arbets- och bruksanvisning
Käyttöohje
Betjenings- og brugsinformation
Μέθοδος και Οδηγίες Χρήσεως

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Method and User Instruction Siloc

The Siloc Bonding System is a method to achieve a chemical - micromechanical bond between veneering materials such as Artglass[®] and Dentacolor[®], or denture base materials and dental alloys.

Applications

The Siloc Bonding System is recommended for use with the following restorations:

Fixed restorations

- Crowns and bridges
- Maryland and Inlay bridges
- Inlays and Onlays
- Implant supported restorations

Removable restorations

- Conus and telescope crowns
- Attachment retained restorations
- Partial dentures
- Metal supported full dentures

Design Features

The Siloc Method was developed to apply a coating onto dental alloys for bond enhancement. The temperature needed to activate the bond coating is generated by indirect heat radiation in the activation chamber of the Siloc unit. A sensor transmits the temperature in the activation chamber as a signal to the electronic control system.

Integrated in the plastic housing is a cooling fan, which avoids overheating of the Siloc unit, even when used in a non-stop mode at high ambient temperatures. There is a heat isolation between the activation chamber and the housing and an isolating flange connects the activation chamber to the unit's frame. The outer housing of the Siloc unit is made from Makrolon (PC) 6557. This material complies with the US safety standard MVSS 302, class V-2 of UL 94 and the FKT regulations. It is self-extinguishing and UV resistant.

Advice

Read the manual first before using the equipment!

Installation

The unit must only be used in rooms with sufficient air ventilation. It must be placed free standing allowing cooling air to circulate. Do not store flammable goods next to the unit. Before connecting the power plug, check the voltage indicated on the ID tag to comply with your supply voltage.

Turn-on

By activating the green main switch, the unit is turned on. The heat-up time is approximately 10 minutes, during which the yellow LED is permanently on, signaling the unit is in "heat-up" mode.

Once the work temperature has been established in the activation chamber, the green LED comes on signaling the "ready" mode. The yellow LED now signals when the heating system is on to maintain the work temperature.

Advice

Switch off the unit only for longer work interruptions or overnight.

Program Selection

Use position 2 of the program dial for all common dental alloys. The program settings 1 and 3 are not used for the Siloc activation.

Correct Framework Design (see graphics page 48)

Always refer to the application chart of the alloy manufacturer. The minimal wall thickness of metal copings is 0.3 - 0.5 mm (see also alloy manufacturer's instructions). When using soft alloys the walls around the cervical margins must be thicker 0.5 mm to ensure sufficient strength.

The metal framework must follow the contours of the tooth in a reduced scale. Depending on the prep design, always replace missing tooth structure with the framework, so that the Artglass veneer can be applied in a uniform layer thickness (cusp support). When total coverage Artglass veneering is applied, make sure the **minimal veneer thickness on posterior occlusions and anterior incisal edges is 1.0 - 1.5 mm!**

Avoid sharp edges and undercuts (exclusive Retention support) on the frameworks. The interproximal surfaces must follow the bow of the tooth shape. Multi-unit bridge frameworks must have strong and non-flexing interproximal connections. Metal reduced pontics must support the Artglass veneer effectively. Fill-in metal reduced areas with Artglass basic. To maintain optimal intraoral maintenance and cleaning, consider periodontic function when designing the framework.

When using retention beads, do not go below the diameter of 0.6 mm.

Alloys

The bond strength on copper-free gold alloys is about 20 % lower than all other gold alloys.

Step by Step Procedure

■ Sandblasting

The veneering surfaces must be blasted with clean, white Al₂O₃ of 110 - 250 micron.

Never reuse the sand! The sandblaster has to be set at a minimal air pressure of 3 bar (42 psi) and the nozzle must be designed to handle the above mentioned grain size.

Advice

Hold the framework during the sandblasting with Kulzer grip or another suitable instrument. A carefully sandblasted alloy surface is an absolute necessary prerequisite for a strong bond, as the metal surface is cleaned and activated. The metal roughening results in a significant increase of active bonding surface and microfine retentions and undercuts are created. It is important to use hoses with the sandblaster that do not get abraded by the sand particles, leaving a contamination layer on the metal surface; always use original replacement hoses made from hard, abrasion resistant plastic material. Check the nozzle regularly for wear and diameter enlargement. A worn nozzle will cause an air pressure drop, which will reduce the acceleration of the sand particles. Consequently, the roughening effect on the alloy surface is reduced, which weakens the bond.

■ Cleaning

The sandblasted surfaces are cleaned only by blowing water and oil free pressured air onto them (30 - 40 psi). Hold the air nozzle about 1 - 2 inches away from the framework.

Advice

If wet-shiny areas develop on the alloy surface during air cleaning, the air contains moisture and cannot be used for cleaning. In case of contamination of the surface, the entire framework has to be sandblasted all over again. Check the oil and water filters of the pressurized air system regularly to maintain absolutely clean and dry air.

■ Application of Siloc pre

Some Siloc pre (white bottle cap and label) is placed in a clean dish; with the white disposable brush, Siloc pre is applied onto the surface as a uniform coating, avoiding puddles and liquid agglomerations.

Siloc pre will wet the rough, sandblasted alloy surface quick and easy. Let air dry (approximately 2 minutes) and place the framework on the Siloc tray.

■ Positioning in the Siloc unit

Place all frameworks (up to 16 units) on the Siloc tray. Position the units outside of the activation chamber. For single crowns and small bridges the use of the positioning pins. Place the units on the cervical margins, not the incisal edges or occlusal tables. By use of the big tweezers, place the Siloc tray in the activation chamber.

Advice

Do not place frameworks directly in the activation chamber, always use the Siloc tray. Use only the original Siloc tray and positioning pins which come with the unit. Do not use other ceramic supports or pins, as the activation procedure will be impaired and as a result the bond fails.

■ Activation

Place Siloc tray in activation chamber, close the door and press "Start" key. Do not keep the door open longer than necessary to place the Siloc tray. If the door is kept open too long, the temperature will drop and the green LED will go off.

In this case, take Siloc tray out of the chamber, close door and let unit heat up again to work temperature. When green LED comes on, place Siloc tray in chamber, close door and start the activation program. The red LED (Program control) indicates the activation program is in process.

After completing the activation, an audio signal comes on. The work pieces have to be removed from the chamber during the time the audio signal is on. When opening the door, the audio signal goes off automatically.

Advice

If the workpieces are not removed during the time the audio signal is on, the signal will go off and the red LED (program control) starts flashing. This indicates that the activation program was overextended and the process has to be repeated all over again, starting with sandblasting. The "Reset" key can be used to switch the unit back to normal function and the flashing LED will go off.

The correct activation can be recognized by a slight discolorization of the metal surface.

■ Program Selection

The program selector is generally set on program 2 for all common dental alloys.

■ Cool Down

Place the framework on the pull-out cooling platform and let cool down to room temperature (approximately 4 minutes).

■ Application of Siloc bond

Place some Siloc bond (red bottle cap and red label) in a clean dish and apply a uniform layer with the red disposable brush. Avoid puddles and liquid agglomerations; let air dry for 5 minutes.

Advice

If Siloc bond is applied on a hot framework, the bond mechanism is destroyed and the procedure has to be repeated all over again, starting with sandblasting. Use only the clean, color coded disposable brushes for the application of Siloc pre and Siloc bond.

Once placed in a dish, the liquids must be used within 5 minutes. Never return liquids from dish, back into the bottles! Immediately after use, close caps of bottles tightly. Use only clean dishes; wipe dishes clean after each use with a paper tissue.

After completion of the activation, the bond surface remains active for approximately 30 minutes.

Never contaminate the surface (finger prints) after sandblasting until the first opaquer layer is placed. Contamination of any kind will act like a separator and causes bond failures. It is, therefore, important to handle the cases in a continuous, uninterrupted procedure, starting with sandblasting all the way to the first opaquer layer.

The luting surfaces of Maryland bridges must be covered with a thin opaquer layer to seal and protect the Siloc layers.

■ Application of Dentacolor Opaquer

Mix powder and liquid 1:1 to a thin consistency. Stir for 30 seconds minimum to achieve a homogeneous distribution of the pigments. Apply a thin coating with a brush and light cure for 90 seconds in the Dentacolor XS or UniXS unit, or BLZ 15 in the Unilux AC. By applying a second and third thin layer and curing each layer separately, a complete masking of the metal surface can be achieved.

■ Application of Artglass Opaque

Place some of the one-component opaquer on the mixing pad, stir it and apply a very thin layer with a short bristle brush. Dab the opaquer into the rough alloy surface; the first layer must not mask the metal yet.

Light cure for 90 seconds (UniXS or Dentacolor XS) and apply a second and, if necessary, a third very thin layer and light cure each layer separately. Each layer should not exceed 75 microns in thickness.

Advice

If the opaquer is not mixed (Dentacolor) and applied properly in very thin layers, the bond is severely impaired! (See user instructions of opaquer for Dentacolor and Artglass too).

Storage Advice

After opening the Siloc pre and bond bottles the first time, they have a shelf life of 6 months. Unopened bottles should be stored in the refrigerator and the shelf life is 2 years. Do not use the materials beyond the expiration date. Once placed in a dish, the liquids must be used within 5 minutes! Always put the caps back on to the bottles immediately after use and tighten firmly. Store away from children's reach. Do not expose to direct sunlight.

Safety Advice

Flammable liquids! Keep containers in a cool and well ventilated place. Do not inhale the fumes! After eye contact, rinse immediately with a large amount of water. Apply antielectrostatic charge measures.

Cleaning and Maintenance

Always disconnect the power plug prior cleaning!

The housing is made of polycarbonate and may be cleaned with a moist cloth and a mild cleaning detergent. Make sure water will not get inside

Attention: Only clean the activation chamber with a mild steel polish; always keep clean and shiny

Never allow wax on any other organic materials into the activation chamber. The pyrolytic fumes will form a deposit on the bonding surfaces and function as a separation medium.

Technical data

Connected load	375 W
Operating voltage (Cycles/Fuses)	100 V ~ 50 - 60 Hz T 5A 120 V ~ 50 - 60 Hz T 5A 127 V ~ 50 - 60 Hz T 5A 230 V ~ 50 - 60 Hz T 2,5A 240 V ~ 50 - 60 Hz T 2,5A
Protection DIN 40050	IP 20
Thermal protection	Class 1
Nominal temperature	250 - 370 °C
The equipment is built in accordance with DIN 12880 standards. If handled in accordance with instructions, health endangering fumes or gases will not be generated.	
Safety rating	Class 1
Purpose	Protection of heat chamber
Extent of protection	In case of malfunction, heat cabinet will not cause danger
Safety provisions according with paragraph 5.5:	Temperature fuses
Safety measures	Special safety measures in relation to intended use
Dimensions	Width approx. 340 mm Height approx. 300 mm Depth approx. 350 mm
Weight	approx. 10 kg

Technical data are subject to change

Liability limitations

The manufacturer is not responsible for other uses the equipment than listed in these instructions; unauthorized changes and tampering with the unit is dangerous and is not covered by the manufacturers warranty and liability. Repairs must only be performed by the manufacturer or authorized service stations by use of original Kulzer repair parts.

Regulatory compliance

The Siloc unit complies with the following safety regulations:

- EC Regulation for low voltage appliances 73/23/EWG
- EC Regulation 89/336/EWG for electro-magnetic compliance

Advice

Keep method and user instruction at an exposed place close the unit!

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